## TRANSPORTING OVER-SIZE PATIENTS



#### **26-B** Heavy Load

- New stretcher with reinforced structure and HEAVY LOAD frame;
- tech. specs: Length (min / max)
   1920mm 1670mm
   Width 550mm
   Loading height
   610/640mm
   Waight 751cm Weight **35kg.**



#### **TSL** Expander Kit

- allows transporting of bariatric patients with SCOOP EXL;
- patented locking system to prevent accidental unlocking
- kit contains:
  2 x extenders (red),
  3 x belt extensions,

# 

Reg. Office & Sales via Zallone, 26 • Pieve di Cento 40066 BO • Italy Phone | +39 051 6860028 Fax | +39 051 6861508

info.it@ferno.com



# NO JOB IS TOO BIG, WITH THE RIGHT EQUIPMENT.

**OFERNO** 





# 

## SUPPORTING OPERATORS IN CHALLENGING HEAVY LOAD RESCUES

Ferno's 26-B Heavy Load and TSL Expander Kit increase the **operator's ability to assist and the comfort of patients**, simplifying the **most complicated and time-sensitive phases of the rescue**.



# Ferno TSL EXPANDER KIT

An robust aluminium accessory that attaches to the Scoop EXL atraumatic trolley. Designed for easy movement and transport of a large patient.

When combined with the Easy Lift system, the Scoop EXL trolley and TSL Expander Kit can also be used with lifting equipment.



**THE KIT CONTAINS:** 2 x red extenders; 3 x belt extensions; 1 x carrying bag.

#### DESIGNED FOR RESCUE WORKERS

**New intuitive handles** for easy control by the operator during the unloading phase.

### REINFORCED STRUCTURE

New lightweight reinforced structure with a **load capacity of 280kg**. Can be reduced in length if needed.

# Ferno 26-B HEAVY LOAD

#### EASY TO USE

Adjustable multiposition backrest and collapsible side rails.

GREATER COMFORT AND SAFETY

Large wheels and new ergonomic brakes designed for use with safety boots.