



EC Declaration of conformity no.: 190523-004 (issue 4)

Product Name:	Vacuum Splint Set Ferno AS 190
Germa Article No:	23005002210
Description:	<p>Vacuum splints are tools for fixation and immobilisation of patients with suspected extremity injuries on arm, leg or hip. The splint consists of an airtight shell with fixation straps, vacuum valve and polystyrene granules all flame retardant.</p> <p>Set of Splints comprising of: 1 x 23005002212, Fore arm splint 1 x 23005002213, Arm splint 1 x 23005002214, Leg splint 1 x 23005002216, Bag for splints 1 x 23005002217, Repair kit 1 x 23272000701, Vacuum Hand Pump</p>
Accessory:	

Identification: All products with serial numbers issued from LOT number: **504040**, date: **2011-05-06** (yyyy-mm-dd).

We the manufacturer hereby declare that this product is in conformity with the requirements in:

- Council directive 93/42/EEC.
Equipment class: I, non-sterile, no measuring function.
- Swedish Medical Devices Act 1993:584
- Regulations for medical devices in LVFS 2003:11

Signed in Kristianstad, Sweden on behalf of AB Germa by:

Position: Managing Director
Name: Björn Holmqvist

Date: 24/05-2019

Sign: 

Postadress	Adress	Telefon	Telefax	Organisationsnummer. F-skattebevis	Bankgiro	Postgiro
Box 2127 S-291 02 KRISTIANSTAD SWEDEN e-mail: info@germa.se	Industrigatan 54-56 KRISTIANSTAD	+46 (0)44-123030	+46 (0)44-103179	556086-9413	481-4125	21 68 22-7