Oximeter







Features

- 1. Display size: 1.5 inch display size, the largest size of fingertip oximeter, Larger numeric font
- 2. Display direction: User-friendly Horizontal Screen Design
- 3. Display assembly process and effect:

The display and the window are completely fitted;

- (1) Display is more transparent and clear.
- (2) larger viewing angle
- 4.PI Parameter: with PI
- 5. Display brightness adjustable: Yes, Three-grades brightness adjustment
- 6.Menu: with Menu, can setup brightness, alarm limits and speaker on/off

Specifications

Model No.	Oximeter: YM101
Dimensions	57mm*31mm*30mm
Weight	Approx.0.3g(exclude 2 AAA batteries)
Display	LED
Display Size	33mm*20mm
Response Time	15 Second
SpO2	
Measurement range	0 to 100%
Display range	0-99%
Resolution	1%
Accuracy	70 to 100%: ±2% 0% to 69%: unspecified.
PR	
Measurement range	25 to 250 bpm
Resolution	1 bpm
Accuracy	±3 bpm
PI	Low Perfusion upto 0.075 %
Power	Battery: LR03*2 Operate current: 25-50mA
Automatic Power Off	11 seconds
Operation temperature	5°C-40°C
Storage temperature	-10°C-40°C
operation	10%-80% in storage
Ctn. / Pcs	48*32*20cm; 100 pcs/CTN. 6.4KGS/CTN