

FERNO Child EasyFix Vacuum Mattress

Designed to ensure effective immobilisation and transportation of paediatrics



**FOR MORE DETAILED
INFORMATION SEE OVER**

Specifications

FERNO – is a world leading innovator and manufacturer of patient transport solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. Ferno's industry-leading **EasyFix Child Vacuum Mattress** typifies our passion for top quality, highly functional design... patient transfer products that don't just do the job, they do it **BETTER**.

Length	Width	Depth	Weight (with straps)	Load Limit
1600 mm	480 mm	45 mm	4.7 kg 5.8 kg (packed in bag with pump & repair kit)	150 kg/23st/330lbs

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk

 **FERNO**[®]

FERNO Child EasyFix Vacuum Mattress

The Ferno Child EasyFix vacuum mattress is the most effective system for the immobilisation of paediatrics.

Incorporating a special body shape design which optimises the vacuum mattress capabilities and minimising the required storage space, the Child EasyFix is designed for the immobilisation of children and infants up to 160 cm and higher, where an effective immobilisation can be achieved.

The Child EasyFix's unique construction and shape delivers unparalleled immobilisation of the upper back, neck and head, reinforced with a base plate along the full length of the cervical spine.

In combination with the vacuum mattress's easy-to-deploy restraint system, and using an X-strap system to support and create the most effective immobilisation, Ferno's Child EasyFix is the optimum system for paediatrics.

Features include:

- Head, neck and upper back shaped area, with integral board for optimum C-spine immobilisation
- Lightweight and easy to store
- X-ray, CT and MR compliant - no need to move the patient to another device for diagnostics
- Suggested patient size from approx 90 cm up to 160cm, but also usable outside this range if an effective immobilisation is achieved
- Kit includes; Vacuum mattress, restraint system, footrest vacuum pump, repair kit and carry bag



Product Code

271504000

Ferno Child EasyFix Vacuum Mattress

271504001

Ferno Child EasyFix Vacuum Mattress Ready2Go Kit (incl. pump, repair kit and bag)

FERNO - Child EasyFix Vacuum Mattress FPC001 04.16

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk



Ferno (UK) Limited, Stubbs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ England.