



5126-EL-IT

Multi-level self-loading stretcher

Code: 5126-EL-IT





5126-EL-IT is the stretcher designed to handle and transport patients in emergency situations outside the hospital. Optimised for use in ambulances, it ensures **high safety standards** thanks to its restraints system and compatibility with Ferno locking systems for rescue vehicles.

MAIN FEATURES

- ▶ **Multi-level.** 5 possible heights facilitate rescue procedures by lightening the load on the back during lifting manoeuvres.
- ▶ **Lightweight and sturdy.** The anodised aluminium alloy structure weighs just 45 kg and can carry up to 250 kg.
- ▶ **Safe.** The rear caster wheels have a safety brake to prevent accidental movements.
- ▶ **Easy to load and unload.** It's equipped with a fall protection system to ensure risk-free manoeuvring during unloading operations.
- ▶ **Can be used as a field stretcher.** The leg locking system keeps the stretcher legs raised and locked, thus allowing it to be used as a field stretcher wherever the terrain does not allow the wheels to be used. The side rails are of the folding type.
- ▶ **Adjustable backrest.** The backrest is adjustable up to a maximum inclination of 80°.
- ▶ **It's designed for performing life-saving manoeuvres** thanks to its rigid ABS tops.
- ▶ **It allows for the transport of incubators**, as it is compatible with the dedicated ITC interface (Incubator Transport Interface).
- ▶ **Easy to clean.** The ABS tops can be removed to make cleaning operations even easier.

AVAILABLE MODELS

DESCRIPTION	PART NO.
Self-loading stretcher, 5 positions, loading height min. 64 cm	5126-EL-IT
Self-loading stretcher, 5 positions, loading height min. 72 cm	5226-EL-IT



COMPATIBLE ACCESSORIES

DESCRIPTION	PART NO.
Four-point restraint	0313915
Two-piece restraint	430-2-N
Safe Transfer Super Comfort mattress	271029000
Pac-Rac (material holder)	0818933
PacRac+	FWEPR+IT
Adapter Kit for PacRac+	25-00214
IV Pole	0087156
Black five-section mattress	25-00200
Slam locking	SLAM-5026/MONDIAL
ITC (incubator fixing interface)	ITC-HL
Neomate (Paediatric transport system)	0314116
Pedimate Plus (Paediatric transport system)	0314115
Patient Shield	0822097
Document holder kit	25-00273

TECHNICAL SPECIFICATIONS

Length	1970 mm
Width	560 mm
Height	(maximum, with backrest up) 1500 mm
Weight	45 kg
Load capacity	250 kg

CERTIFICATION



The product meets the general safety and performance requirements of **Regulation (EU) 2017/745 (MDR)** concerning medical devices.

Compliant with **EN1865** and compatible with the SLAM locking system for ambulance transport (**EN 1789**-certified system).