

WIZLOC JUNGLE

Multi-sized Pediatric Cervical Collar

Code Number 21-0140-020



WIZLOC JUNGLE



Ferno's WizLoc JUNGLE collar is a multisized cervical immobilizer specifically studied for traumatized patients in the age bracket including 2 to 10 years. The pediatric Wizloc Jungle collar allows you to manage all types of injuries resulting from head trauma.

FEATURES

- *Mono-valve use: easy and quick positioning.*
- *Bi-valve use: can be adapted to patients with particular needs.*
- *Can be placed into position either from the left or right side of the patient according to the patient or operator needs.*
- *Micrometric front adjustment: major adaptability to the patient and neutral spine alignment.*
- *Rear section adjustment: can be used on an ample range of pediatric patients without occipital compression.*
- *Adjustment directly on the patient with collar in place.*
- *Rear section can easily be shaped and adapted to any form of the occipital area.*
- *Hypoallergenic lining, closed cell, to avoid absorption of biological liquids.*
- *Totally x-ray transparent and MRI compatible.*
- *Can be totally decontaminated (up to 10 times).*
- *Numbered decontamination stickers to affix directly on the collar (provided).*
- *Behavior instructions for parents translated into the most common languages spoken.*
- *Ample range of use: pediatric patients from 2 to 10 years of age.*
- *Training course on pediatric trauma management from the Ferno Academy.*

TECHNICAL SPECIFICATIONS

Max Length	open: 450 mm
Height: front	160 mm
Height: back	150 mm
Weight	180 g
Range of application	2-10 years old

Code

21-0140-020

Description

WizLoc JUNGLE, multisized pediatric cervical collar, comes with adhesive numbers to keep track of how many times it has been decontaminated and with an album and stickers for children.